

# VIKING-EC200

## Laboratory Conductivity Meter · EC / TDS / Temperature

Model: VIKING-EC200  
 Product number / SKU: PRO2541  
 Product family: Conductivity meters  
 Format: A4 landscape · EN

**VIKING-EC200** is a laboratory conductivity meter for routine measurement of **electrical conductivity (EC)**, **total dissolved solids (TDS)** and temperature in water quality, quality control, process and teaching laboratory workflows.

The instrument combines a **6.0" high-resolution LCD with backlight**, Auto Read and Continuous measurement modes, automatic/manual temperature compensation, 1-point calibration with standard recognition and GLP-style result storage.

EC / TDS / Temp

6.0" LCD

±1.0% FS

1-point calibration

ATC / MTC

50 results / parameter

IP54

### Conductivity range

0.00  $\mu\text{S}/\text{cm}$  to 200  $\text{mS}/\text{cm}$  with range-dependent resolution down to 0.01  $\mu\text{S}/\text{cm}$ .

### Standard kit

PRO6154 EC electrode, 1413  $\mu\text{S}/\text{cm}$  standard solution and EHO5 electrode stand.

### Laboratory workflow

Auto Read, Continuous reading, Auto-Hold prompts and temperature compensation.

### Online reference

Open product page · Electrochemical instruments



## Key features and supplied configuration

**Routine EC / TDS measurement:** VIKING-EC200 is intended for laboratories that need a simple and reliable benchtop conductivity meter with temperature compensation, guided stability prompts and compact GLP-style data storage.

### Key instrument features

- 6.0" high-resolution LCD with backlight.
- Multi-reading modes: Auto Read and Continuous.
- Automatic and manual temperature compensation: ATC / MTC.
- Auto-Hold function with on-screen prompts: Reading / Stable / Locked.
- GLP-style data storage: 50 results for each parameter.
- 1-point calibration with automatic standard recognition.
- Settable cell constant and factory reset to default settings.
- IP54 protection for routine laboratory environments.
- Auto shutdown options: 300 / 600 / 1200 / 1800 / 3600 s, or Off.

### Measurement scope

- **Conductivity:** 0.00  $\mu\text{S}/\text{cm}$  to 200 mS/cm.
- **TDS:** 0.00 ppm to 100 ppt.
- **Temperature:** -5.0 to 110 °C.
- **Accuracy:**  $\pm 1.0\%$  FS for conductivity and TDS.
- **Reference temperature:** 25 °C.
- **Inputs:** conductivity probe with temperature sensor, 5-pin aviation connector.

### Included items

Conductivity electrode **PRO6154**

1413  $\mu\text{S}/\text{cm}$  EC standard solution, **50 mL**

EHO5 electrode stand **PRO6110**

Main instrument **VIKING-EC200**

AC adapter, 100–240 V AC input

DC 9 V output power supply

### Routine water quality

Measurement of conductivity, dissolved solids and temperature in basic laboratory water analysis.

### QC laboratories

Stable benchtop operation for incoming water, process water and quality control checkpoints.

### Teaching labs

Clear display, simple calibration and limited parameter set for training and repeatable routines.

## Technical specifications

### Conductivity and TDS

<b>Model</b>	VIKING-EC200
<b>Measured parameters</b>	EC / TDS / Temperature
<b>Conductivity range</b>	0.00 $\mu$ S/cm to 200 mS/cm
<b>Conductivity resolution</b>	0.01 $\mu$ S/cm minimum, varies by range
<b>Conductivity accuracy</b>	$\pm$ 1.0% FS
<b>Reference temperature</b>	25 °C
<b>Calibration points</b>	1 point
<b>Standard recognition</b>	84 $\mu$ S/cm, 1413 $\mu$ S/cm, 12.88 mS/cm
<b>TDS range</b>	0.00 ppm to 100 ppt
<b>TDS resolution</b>	0.01 ppm minimum, varies by range
<b>TDS accuracy</b>	$\pm$ 1.0% FS

### Inputs, protection and power

<b>Inputs</b>	Conductivity with temperature probe; 5-pin aviation connector
<b>IP rating</b>	IP54
<b>Power</b>	AC adapter, 100–240 V AC input; DC 9 V output
<b>Dimensions</b>	242 × 195 × 68 mm
<b>Net weight</b>	900 g / 1.98 lb

### Temperature, reading modes and memory

<b>Temperature range</b>	–5.0 to 110 °C
<b>Temperature units</b>	°C
<b>Temperature resolution</b>	0.1 °C
<b>Temperature accuracy</b>	$\pm$ 0.2 °C
<b>Reading modes</b>	Auto Read, Continuous
<b>On-screen prompts</b>	Reading / Stable / Locked
<b>Temperature compensation</b>	ATC, MTC
<b>Data storage</b>	50 results for each parameter
<b>Backlight</b>	Yes
<b>Auto shutdown</b>	300 / 600 / 1200 / 1800 / 3600 s; Off

#### Calibration workflow

- Prepare the conductivity electrode and calibration standard.
- Use one-point calibration with recognized standards.
- Verify stable reading prompt before recording results.
- Use ATC or MTC according to the laboratory method.
- Store selected results in the instrument memory for GLP-style traceability.

## Laboratory use, selection notes and comparison context

### When VIKING-EC200 is a suitable choice

- Routine conductivity and TDS checks where a compact benchtop meter is preferred.
- Laboratories that need a simple display, basic calibration and temperature compensation.
- Teaching laboratories and QC rooms where workflows should remain easy to repeat.
- Water quality screening, process water monitoring and incoming water checks.

### When to consider a higher model

- If resistivity, salinity or wider conductivity range is required.
- If multi-point calibration is important for the method.
- If larger GLP memory, advanced user management or touchscreen operation is required.
- If USB/RS-232 connectivity is needed for data export or printer connection.

### Important measurement considerations

Conductivity measurement depends on the electrode cell constant, temperature compensation settings, calibration standard condition, sample matrix and cleanliness of the probe. For routine water analysis, the 1413  $\mu\text{S}/\text{cm}$  standard is a practical calibration point, but the final method should define the expected conductivity range and acceptance criteria.

The supplied PRO6154 EC electrode and EHO5 stand make the instrument suitable for stable benchtop use. For more complex water chemistry, ultra-low conductivity, high-salinity samples or extensive documentation requirements, the configuration should be reviewed before quotation.

### Method fit

Best for EC, TDS and temperature routines rather than advanced multiparameter analysis.

### Documentation

50 results per parameter support simple GLP-style result retention.

### Bench layout

A compact instrument footprint and electrode stand support clean benchtop workflows.

## Product links, video placeholder and technical notice

### Configure VIKING-EC200 for routine EC and TDS laboratory workflows

COLO.Science can support selection of conductivity meters, electrodes, calibration standards and accessories for water quality, QC, teaching and routine laboratory measurement.

[Open product page](#)[Benchtop conductivity meters](#)[Electrochemical instruments](#)[Request quotation](#)[PDF version](#)

#### COLO.Science

Laboratory equipment, technical specifications and product configuration.

sales@colo.si  
+386 64 222 724  
<https://colo.si/contacts/>

#### Recommended document links

- **Product page:** VIKING-EC200 product page
- **Category:** Benchtop conductivity meters
- **Electrochemistry page:** Electrochemical instruments
- **HTML path:** /wp-content/images/TechSpec/WaterQualityP/VIKING-EC200.html
- **PDF path:** /wp-content/images/TechSpec/WaterQualityP/VIKING-EC200.pdf

#### Video presentation · in preparation

A VIKING-EC200 product video is planned. It will present the display, calibration workflow, EC/TDS measurement routine, electrode connection and standard delivery configuration.

**Temporary video placeholder:** #viking-ec200-video

Replace this anchor with the final YouTube / product video URL when available.

**Technical disclaimer:** This COLO.Science technical specification is prepared for product orientation, quotation preparation and preliminary comparison only. Technical characteristics, supplied accessories and configuration options may vary depending on the confirmed offer, manufacturer documentation and selected electrode configuration. Only the official quotation, order confirmation and manufacturer-approved specification should be treated as definitive for procurement, tender or contractual use.