



CENT600C Textile Balance with Cutter

External Calibration

Balance Specification

- Balance Capacity 600g
- Resolution 0,01g
- A High quality load cell
- Full range tare
- Overload protection
- Under hook
- Antivibration Filter

GSM Cutter specification

- Apply for textile fabric, carpet, film, foam, paper and board
- Cutter area 100cm²
- Cutting Blades: Reversible High quality steel sharp and durable
- Sample diameter: 113cm
- Standard: GB/T4802.1 and GB/T4669
- Cutting depth: up to 5mm
- Application: Textile fabric, carpet, film, foam, paper and board

Download
Catalogue



CENT600c

High Capacity Precision Balance



External Calibration

External Linear Calibration




High Precision Load Cell

High measuring range up to 30kg is provided with high capacity load cell,



Software for data collection

Optional software for communication with the scale provides the ability to collect data, control the measurement process ...



GSM Cutter 100cm²

600g/0,01g

External
Calibration

Big Pan Size

143*192mm

High Precision
Load Cell

Piece Counting