




## LACE18 Laboratory High Speed Centrifuge



LED Display



Neutral  
Temperature

Download  
Catalogue



Price - €



### 18500rpm/23990xg

Maximum speed and rcf meet all demand for research grade centrifuge



### LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work



### Lead Lock

Simple Lead Lock with one button motorized



### Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



### Safety Function

Automatic Electric Lead Lock with alarm with safety function, disbalance detection, over speed protection

Item Name	LACE18
Max. Speed rpm	18500rpm/23990xg
Max. Capacity	4x800ml
Timer	0-99min
Speed Accuracy	±20rpm
Display	LED – Display: Speed, Time, RCF
Noise	≤65dBA
RPM/RCF conversion	Yes
End process sound alarm	Yes
Electric lid lock	Yes
Safety	Mechanical lid lock
Power supply	AC 220V/50Hz
Dimension	685x500x385
Net Weight (kg)	70

LACE18 Rotors				
Rotor NO	Type	Max Speed/RCF	Max. Volume	Test tubes
PRO	Swing Out Rotor	4000rpm/3450rcf	4x800ml	800ml
PRO2830	Microplate Rotor	4000rpm/2840rcf	4x4x96well	
PRO	Microplate Rotor	4000rpm/2490rcf	2x4x96well	
PRO	Microplate Rotor	4000rpm/2300rcf	2x3x48well	
PRO	Swing Out Rotor	4000rpm/2840rcf	4x30x5ml vacuum tube	5ml vacuum tube
PRO		4000rpm/3140rcf	4x30x7ml vacuum tube	7ml vacuum tube
PRO		4000rpm/3140rcf	4x18x10ml vacuum tube	10ml vacuum tube
PRO	Swing Out Rotor	4000rpm/3830rcf	Fat bottle	
PRO	Angle Rotor	16000rpm/15760rcf	4x8 PCR	
PRO	Angle Rotor	15000rpm/21420rcf	6x8 PCR	
PRO	Angle Rotor	16000rpm/17480rcf	8x8 PCR	
PRO	Angle Rotor	15000rpm/22930rcf	12x8 PCR	
PRO	Angle Rotor	15000rpm/22920rcf	40x0.5ml	0.5ml
PRO	Angle Rotor	21000rpm/30910rcf	12x1.5/2ml	1.5/2ml
PRO	Angle Rotor	16000rpm/23440rcf	24x1.5/2ml	1.5/2ml
PRO2841	Angle Rotor	14000rpm/20800rcf	30x1.5/2ml	1.5/2ml
PRO	Angle Rotor	13000rpm/17930rcf	48x1.5/2ml	1.5/2ml
PRO	Angle Rotor	16000rpm/22020rcf	16X5ml	5ml
PRO	Angle Rotor	16000rpm/21500rcf	6x10ml	10ml
PRO	Angle Rotor	15000rpm/22680rcf	12x10ml	10ml
PRO	Angle Rotor	13000rpm/19490rcf	16X10ml	10ml
PRO	Angle Rotor	13000rpm/17790rcf	8x15ml	15ml
PRO	Angle Rotor	11000rpm/14330rcf	12x15ml	15ml
PRO	Angle Rotor	5000rpm/3500rcf	24x15ml	15ml
PRO4471	Angle Rotor	5000rpm/3830rcf	30x15ml	15ml
PRO	Angle Rotor	14000rpm/19060rcf	6x30ml	30ml
PRO	Angle Rotor	13000rpm/18840rcf	6x50ml conical	50ml conical

Please note

If this Centrifuge does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

## COLO LabExperts

Laboratory Equipment production and Distribution

www.colo.si

Presernov Trg 9

8000 Novo Mesto

Slovenia

tel. +386 590 48 880 E-mail: office@colo.si



www.colo.si

© 2021 COLO Distributors Price List V1

LACE18 Rotors				
Rotor NO	Type	Max Speed/RCF	Max. Volume	Test tubes
PRO	Angle Rotor	13000rpm/18730rcf	6x50ml round	50ml round
PRO	Angle Rotor	5000rpm/3860rcf	12x50ml	50ml
PRO	Angle Rotor	4000rpm/2970rcf	24x50ml	50ml
PRO	Angle Rotor	13000rpm/18940rcf	4x85ml	85ml
PRO	Angle Rotor	12000rpm/14850rcf	4x100ml	100ml
PRO	Angle Rotor	10000rpm/11380rcf	6x100ml	100ml
PRO	Angle Rotor	4000rpm/2970rcf	12x100ml	100ml

Please note

If this Centrifuge does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

## COLO LabExperts

Laboratory Equipment production and Distribution

www.colo.si

Presernov Trg 9

8000 Novo Mesto

Slovenia

tel. +386 590 48 880 E-mail: office@colo.si



[www.colo.si](http://www.colo.si)