



LED Display



Neutral  
Temperature

Download  
Catalogue



## Max. Capacity 4x250ml

Swing out Rotor



## 5000rpm/4807xg

Maximum speed and rcf meet all demand for clinical centrifuge



## TOUCH Screen Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



## Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



## Safety Function

Electric Lead Lock with safety function, disbalance detection, over heating of centrifuge protection

# LACE32T Low Speed Centrifuge

## Angle/Swing out Rotor

- 30569 - Swing out, 4x250ml - 4000rpm/2880xg
- 31376 - Swing out, 24x15ml - 4000rpm/2950xg
- 31377 - Swing out, 32x10ml - 4000rpm/2880xg
- 31378 - Swing out, 48x5 ml - 4000rpm/2480xg
- 31379 - Swing out, 48x7ml - 4000rpm/2760xg
- 5567 – Swing out, 8x100ml - 4000rpm/3077xg
- 31380 - Swing Micro plate, 2x2x96well -
- 30641 - Angle, 12x50 - 5000rpm/3860xg
- 30639 - Angle, 24x15ml - 5000rpm/3500xg

Item Name	LACE32T
Max. Rotor capacity	4X250ml
Max. Speed rpm/xg	5000rpm/4807xg
Power	450W
Ambient temperature	0-40°C
Noise	≤55dB (A)
Timer	0-99h59nin / continuous / short spin button (option)
Electric safety lid lock	Mechatronic motor door lock
RPM/RCF conversion	Yes
Acceleration/Brake	10 programming step of acceleration and deceleration
Power supply	AC 220V/50Hz
Weight kg NV/GW	37/41