



## LACE-5A Low Speed High Capacity Universal Centrifuge



LED Display



Neutral  
Temperature

Download  
Catalogue



Angle/Swing out Rotor

LACE32 is an economical universal centrifuge suitable for large clinical laboratories.

30690 - Swing out, 72x10ml - 4000rpm/3140xg

31377 - Swing out, 32x15 ml - 4000rpm/3200xg

31379 - Swing out, 88x7/5ml - 4000rpm/2760xg

31461 - Swing Micro plate, 2x4x96well 4000rpm/2490xg

30641 - Angle, 12x50 - 5000rpm/3860xg

30642 - Angle, 24x50 5000rpm/3860xg

30614 - Angle, 6x100ml - 5000rpm/3500xg

30614a - Angle 12x100 - 5000rpm/3500xg



### Max. Capacity 4x300ml

Swing out Rotor



### 5000rpm/4830xg

Maximum speed and rcf meet all demand for clinical centrifuge



### Swing out/Angle rotors

Adapters available for this rotor



### LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



### Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



### Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

Item Name	LACE-5A
Max. Rotor capacity	4x300ml
Max. Speed rpm	5000rpm
Noise	<58dB (A)
Timer	0-99 min
Safety lid lock	Mechanic door lock
Power supply	AC 220V/50Hz
Product dimension	620 x 450 x 355mm
Weight kg NV/GW	52kg/60kg

[www.colo.si](http://www.colo.si)

© 2020 COLO Catalogue V.1.0 Rights Reserved.