



LACE-5AH Low Speed High Capacity Universal Centrifuge



LED Display



Neutral Temperature

Download Catalogue



Angle/Swing out Rotor

- Brush-less frequency motor drive, maintenance is free.
- Table top low-speed multipurpose centrifuge.
- Melt housing, 304# stainless steel centrifugal chamber, safe and reliable.
- Microprocessor control and digital display run parameters, like speed, RCF, timer.
- Brushless motor
- Lid lock safety
- Acceleration and breaking adjustment
- Steel rotor housing and metal centrifuge body
- Overheating and over-weight protection
- Conversion of RPM and RCF
- Display shows rotation speed, time and RCF

LAC5AH is High Capacity, economical universal centrifuge suitable for large clinical laboratories.



Max. Capacity 4x800ml

Swing out Rotor



5000rpm/3860xg

Maximum speed and rcf meet all demand for clinical centrifuge



Swing out/Angle rotors

Adapters available for this rotor



LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

- 30696 Swing out Rotor 4 x 500ml 4000rpm/3380xg
- 30671 Swing out Rotor 4 x 800ml 4000rpm/3450xg
- 31494 4x4x96 Micro plate Rotor