



LACE-5RS Refrigerated Low Speed High Capacity Centrifuge

Angle/Swing out Rotor



LED Display



Refrigerated
-20°C to 40°C

Download
Catalogue



- 30569 Swing out Rotor 4x250ml 4000rpm/2880xg
- 3277 Swing our Rotor 32x10ml 4000rpm/2880xg
- 31376 Swing out Rotor 24x15ml 4000rpm/2950xg
- 31378 Swing out Rotor 48x5ml 4000rpm/2480xg
- 31379 Swing out Rotor 48x7ml 4000rpm/2760xg
- 31380 Micro plate Rotor 2x2x96Well 4000rpm/2490xg
- 30614 Angle Rotor 6x100ml 5000rpm/3130xg
- 30641 Angle Rotor 12x50ml 5000rpm/3860xg



Max. Capacity 4x250ml

Swing out Rotor



5000rpm/4730xg

Maximum speed and rcf meet all demand for clinical centrifuge



Swing out/Angle rotors

Adapters available for this rotor



LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

Item Name	LACE-5RS
Max. Speed (Swing out rotor)	5000rpm/4730xg
Max. Capacity (Swing out rotor)	4 x 300ml
Noise	<58dBA
Timer	99H, 59min
Rotor type	Swing out and Angle Rotor
Power supply	AC 220V/50Hz
Product dimension	650 × 665 × 400mm
Weight kg NV/GW	118kg/130kg