



MICE12S Small Laboratory Centrifuge 5000rpm/3600xg

Little Giant

 Analog Display

 Neutral Temperature

[Download Catalogue](#) 

Laboratory centrifuge mice12s with the 12t-place fixed rotor, is an economical solution for smaller laboratory work. Small size and weight allow for easy transportation from one place to another. MICE12S Centrifuge is ideal as an additional instrument in a laboratory, serving either as a reserve or for sample preparation.

MICE12S Small Laboratory Centrifuge light weight, small size, proper speed for blood , urine PRP separation...



5000rpm/3600xg

Maximum speed and rcf meet all demantd for clinical centrifuge



Angle Rotor 12x20ml

Combined with adapter this type, or rotor can work with all type of centrifuges from 2ml to 20 ml,



ANALOG Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Free maintenance Motor



Lead Lock

Simple Lead Lock with one button

PRODUCT CODE	LACE12S
PRICE	
Speed	5000rpm/3600xg
Capacity	12x20ml
Timer	0-99min
Noise	<45dBA
Display	ANALOG – Display
Power supply	220V (±10%) 50-60Hz
Safety	Mechanical lid lock
Weight NW/GW kg	7/8,5kg
Product Dimension	260×265×310mm