





MICE8 Small Laboratory Centrifuge 4000rpm/2800xg


Little Giant


-  LED Display
-  Neutral Temperature


Download Catalogue 


Laboratory centrifuge MICE8 with the eight place fixed rotor, is an economical solution for smaller laboratory work. Small size and weight allow for easy transportation from one place to another. MICE8 Centrifuge is ideal as an additional instrument in a laboratory, serving either as a reserve or for sample preparation.


MICE8 Small Laboratory Centrifuge light weight, small size, proper speed for blood , urine PRP separation...

 **4000rpm/2800xg**
Maximum speed and rcf meet all demand for clinical centrifuge

 **Angle Rotor 8x15ml**
Combined with adapter this type of rotor can work with all type of centrifuges from 2ml to 15 ml,

 **Digital Display**
Speed/RCF, Time, Parameters can be set during the centrifuge work.

 **Brushless Motor**
COLO brushless motor is maintenance free, no carbon dust pollution

 **Lead Lock**
Electromagnetic lead lock

| | |
|-------------------|---|
| PRODUCT CODE | MICE8 |
| PRICE | |
| Speed | 4000rpm/2800xg |
| Capacity | 8x15ml |
| Timer | 0-99min |
| Noise | <45dBA |
| Display | LED – Display: Speed, RCF, Time, RPM/RCF conversion |
| Power supply | 220V (±10%) 50-60Hz |
| Safety | Electromagnetic lid lock |
| Weight NW/GW kg | 7/8,5kg |
| Product Dimension | 260x265x310mm |