



MINICEN7T Minicentrifuge with multi rotor

MINICEN7T

Download
Catalogue



Multi rotor




Small size - large capacity




High centrifugal power



Ergonomic design



Type	MINICEN7T	
Control mode	PI direct drive – until press lid switch is on	
Timing	No	
Display	No	
Max relative centrifuge force	3286xg	
Max rotation speed	7000rpm	
Max speed acceleration/reduction	≤12s / ≤15s	
Noise	≤47db(A)	
Input power		
Power supply	AC100-240V / DC24V	
Ambient conditions	5-40 °C / 80%RH	
Rotor specification		
Rotor type	Picture	Rotor capacity
Multi rotor included		0.2/0.5/1.5/2ml micro tube*12 + 0.2ml 8 PCR stripe*4 (32*0.2ml)

Please note

If this Mini centrifuge does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will produce technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts



Laboratory Equipment production and Distribution

www.colo.si

Polje ob Sotli 4

SI-3255 Slovenia

tel. +386 590 48 880 E-mail: office@colo.si