



DEN100 Densitometer

Download
Catalogue



- Measurement Range^{**}: Typically from 0.00 to 15.000McFarland units³.
- ^{**}Measurement Time^{**}: Approximately 1 second³.
- ^{**}Precision^{**}: $\pm 3\%$ ³.
- ^{**}Standard Deviation^{**}: At 3.00 McFarland units, the standard deviation is ± 0.1 McF³
- Volume sample: >2ml

- Tube Holder: Supports 18 mm glass tubes with included adapter 16mm
- Display: LCD/0.01 readings in McFarland units.
- Working temperature range +4 to 40°C
- Power Supply: Can be powered by an external power supply or three AA batteries
- Set of standards: included
- Dimension: 75*115*165
- Weight: 0.7kg
- LED light source
- Wavelength: 565nm

COLO LabExperts
Polje ob Sotli 4
SI-3255 Slovenia
www.colo.si
Email: sales@colo.si