

EPVTCT4 Vertical Electrophoresis Tetra Cell Tank

Product Details

Manufacturer: Colo Lab Experts

COLO CODE: EPVTCT4



Technical Specification

- Installation for only 15 second, without knob fast and convenient
- Auto switch off when the lid is removed
- Buffer provide cooling effects and keeps pH value stable during the whole experiment process
- Background color provide easy observation of sample adding and electrophoresis process

Technical specification:

- Glass plate (W*L) : 100*100 mm
- Gel size (W*L): 83*73 mm
- Total buffer volume: 1000 ml (for two gels 700-800 ml)
- Gel thickness 1,0 and 1,5 mm
- Number of gels: 1-4
- Number of samples: 11-60
- Volume (L*W*H) : 220*170*120 mm
- Weight: 2,5 kg

Packing:

- Body tank with electrode 1 PC
- Gel casting base 4 PC
- Notched glass plates 6 PC
- 1,0 mm bonded glass spacer 4 PC
- 1.5 mm bonded glass spacer 4 PC
- 1,0 thickness 11 wells comb 4 PC
- 1,5 thickness 11 wells comb 4PC
- 1,0 thickness 15 wells comb 4 PC
- 1,5 thickness 16 wells comb 4PC
- Gel Shovel 4 PC
- Clamp version 2 PC
- Dummy plate 3 PC
- Electrical cables set 1 PC

Optional Accessories