

DBI10 Dry Bath Incubator RT +5 to 100C Main Body

Product Details

Product Number

PRO4538



Technical Specification:

Product model	DBI10	DC10
Temp setting range	0□~100□	-10□~100□
Temp control range	RT+5□~100□	RT- 30□~100□
Time range	1min ~ 99h59min/∞	
Temp accuracy	≤±0.3□	
Display accuracy	0.1□	
Temp uniformity	≤±0.3□	
Heating time	≤12min(from 25□ to 100□)	≤12min(from 25□ to 100□)
Cooling time	-- --	≤10min(from 100□ to 25□)
Cooling time	-- --	≤25min(from RT to RT-30□)
Storage at 4□	-- --	Yes
Multi points running	Yes (Max 5 points)	
Multi cycle	Yes (Max 99 times)	
Auto running	Yes	
Auto heating	Yes	
Auto resume to run	Yes	
Capacity	See optional block	
Power	150W	
Power input	AC110 or AC220V / 50-60HZ	AC100~240V, 50/60Hz
Fuse	250V 3A Φ5×20	
Size	260x195x150mm	
Weight	2.5Kg	4.5Kg

Accessories: