



BOTFIX Bottle top Dispensers

Download
Catalogue



According to ISO 8655 standards



Optional Calibration 17025 ISO



Chemical resistant



High accuracy



Fixed volume

- Developed innovative equipment to ensure the easy and safe handling of corrosive reagents

Type	BOTFIX2.5	BOTFIX5	BOTFIX10
Volume ml	2.5	5	10
Accuracy ±mL	0.0125	0.025	0.05
CV ±mL	0.005	0.01	0.01
Measurement error	According to ISO 8655-5		
Autoclaving	fully autoclavable at 15 psi and 121°C for a duration of 15-20 minutes		

COLO Laboratory Equipment
Polje ob Sotli 4
SI-3255 Slovenia
Tel.+386 64 222 724
Email: sales@colo.si
www.colo.si