



LHDM900D Laboratory Emulsification Homogenization Digital

With standard tools head,

Price - €

Download
Catalogue



The LHDM900D Homogenizer is an ideal laboratory device for higher viscosity samples. Digital control, selection of rotors and stators according to your needs, will provide your work with the necessary tools.



200rpm-11000rpm

Digital Speed Regulation



Maximum Viscosity Mixing 12000Mpas

Ideal for pharmaceutical laboratory viscosity grade



Maximum mixing volume 40L

High laboratory mixing volume ..



Different shape of tools

Depand of the viscosity LHDM900 use diffeent shape and diameter tools



Digital control

To ptovide standard in your work

Item Name	LHDM900D Homogenizer
Max. Speed rpm	200-11000
Max. Volume capacity	40 L
Max Viscosity	12000Mpas
Torque	34.1N.cm
Stator Size	50mm/20mm/5mm
Continually Working	3-10min
Dimension (L*W*H)	650*250*720mm
Power	500W
Power supply	AC 220V ±10%-50/60HzHz
Weight kg GW	18
Packing size	540mm*520mm*200mm

Please note

If this Homogenizer does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts



Laboratory Equipment production and Distribution

www.colo.si

Presernov Trg 9

8000 Novo Mesto

Slovenia

tel. +386 590 48 880 E-mail: office@colo.si

www.colo.si

© 2020 COLO Technical Specification