

OVST – Overhead Stirrer Digital/Analogue



Download
Catalogue



Designing for continuous operation



Designed to save space in your lab



Constant speed under changing loads




Constant torque over all range



Analogue/Digital Display



Stirrer type	OVSTA20	OVSTD20	OVSTD40
Price			
Stirring control	Analog	Digital	Digital
Maximum stirring volume L	20L	20L	40L
Speed RPM	100-2500	100-2500	100-2500
Timer	X	Yes	Yes
Driving mode	Brushless DC motor, maintenance free		
External size mm	115×95×240		
Plate size mm	200*300		
Maximum plate load	50kg		
Voltage	110-230V 50-60Hz		
NW/GW kg	5/6	6/7	7/8



Please note
 If this Stirrer does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts
 Laboratory Equipment production and Distribution
 www.colo.si
 Presernov Trg 9, 8000 Novo Mesto, Slovenia
 tel. +386 590 48 880 E-mail: office@colo.si

