



VOMI25 Vortex Mixer

Price - €

Download
Catalogue



Analogue Speed Regulation

0-3000 rpm



Continuous /Touch rotating operation



DC brushless motor

Maintenance free



Orbital shaking motion



Compact design

- The Vortex Mixer provides smooth instant vortex of even the largest sample sizes.
- Available for different turn plate and laboratory containers.
- Functions of vortex mixing - orbital mixing
- Orbital diameter 4 mm.
- Rotate speed by electronic adjustable
- Touch or continuous operation mode
- Multi-function Vortex Mixer
- Includes a standard top head for mixing a single test, the possibility of replacing the adapter
- Speed electrical controlled 0 - 3000rpm,
- Speed display: scale
- Motor type: Shaded pole motor
- Pop-off Cup Head for test tubes up to 50ml (30mm)
- Shaking Motion: Orbital
- Control: Analogue Phase Controller
- Dimensions (WxHxD) (mm): 127x130x160
- Permissible ambient temperature: 5-40°C
- Permissible relative humidity: 80% RH
- IP 21 level protection
- Body: robust aluminium-cast construction
- Steel base with specially designed vacuum suction
- Silicone feet provides stability and dampens vibration
- Chemical resistant
- Weight: 3,5 kg
- Power: AC 200-240 V, 50 / 60 Hz - 60W