

UCDCH3 Ultrasonic Cleaner digital control with heating 3L

Product Details

Product Number

PRO2471



Technical Specification:

- Large-screen illuminated LED display for time and temperature
- Volume: 3 l, with digital timer and heating
- Parameter selection menu on the display
- Power of 40% -100% adjustable
- Digital ultrasonic cleaning time setting
- Stainless steel
- Heaters made of Kanthal
- Working time countdown display
- Operating parameters of power and memory
- Adjustable heating up to 80C
- WiFi temperature monitoring and control
- Timer range: 1-99 minutes
- Output frequency: up to 145 kHz
- With a special stainless steel basket, sound reduction lid
- Instruments and sound reduction lid inner and outer shells made of stainless steel
- Housing made of stainless steel
- Operating procedures of the instrument using SCM software
- Tank dimensions (L x W x H): 240 x 140 x 100 mm
- External dimension (L x W x H): 270 x 170 x 240 mm
- Power supply: 230 V / 50 Hz
- Power W: 180
- Heating power W: 400
- Drain: Yes

Accessories: