

24 port SPE Glass Block Vacuum Manifold,

Product Details

Product Number

PRO5905



Technical Specification:

1. **Glass Block Design:** The manifold features a glass block lid that allows visible inspection of the extraction process.
2. **Vacuum Control Valve:** You can optimize sample flow rates during extraction using the vacuum control valve.
3. **Test Tube Racks:** These manifolds support 12-13mm or 16mm disposable test tubes, allowing flexibility in sample collection.
4. **Vacuum Gauge and Bleed Valve:** Easily monitor the vacuum process with the built-in gauge and bleed valve.
5. **Durable and Chemically Resistant:** The Glass Block Manifolds are chemically resistant and designed for long-lasting use.

For specific configurations:

- **24-port Vacuum Manifold:** Ideal for processing multiple samples simultaneously.
- **16-port Vacuum Manifold.**

Accessories: