

# PORTpH50 Portable pH/Temp. mV/ Meter 3 in 1 electrode

## Product Details

Product Number

PRO6078



## Technical Specification:

- Measuring Range: pH: (2.00~20.00 pH)
- Range: 2.00 - 20.00
- Resolution: pH: 2,000.-19,999
- mV:  $\pm 1200$   
 $\pm 2500$
- Temperature range: 5.0 to +105 °C
- Accuracy: pH:  $\pm 0.1$ pH  
mV:  $\pm 0.3$  mV  
Temperature:  $\pm 0.1$ °C
- Measuring: Manual or timing controlled
- Display: LCD graphical with backlight
- Automatic/Manual temperature compensation ensures accurate results
- Sentix 3 in 1 pH electrode with integrated temperature sensor
- 1-3 points calibration with Standard Recognition.
- Data capacity of up to 200 sets for each parameter
- Auto-power off feature effectively extends the battery service life.
- IP65 waterproof.
- Carrying case with attachment for accessories
- pH buffers 4.01, 7.01, 10.01
- Continuous working mode more than 850 hour without backlight and 120 h with backlight
- USB cable
- Holder
- Measuring beaker
- 4xAA ACCU battery

## Accessories: